



R82<sup>®</sup>

# R82 Crocodile

Supporting the motivation to move

# Supporting the motivation to move

The R82 Crocodile is a walking aid for children and youngsters with disabilities who need supplementary assistance with their daily walking. The R82 Crocodile is a highly adaptable walking aid for children and youngsters, that encourages movement and more independence, while increasing participation and engagement. It is suitable for children with a body height from approximately 60 cm.

The Crocodile helps the child into a more upright posture. It is placed behind the child and its lightweight design makes it easy to move and use. As much energy as possible is used walking rather than moving the frame.

The multi adjustable handles can be positioned and adjusted for both height and width. Positioned close to the body, they can give better support and stability.

The Crocodile offers a wide range of accessories. The accessories can provide the right amount of support for the user to feel safe and independent.



## Experience the Crocodile in a 360° view



# Features

Individual product solutions can be achieved by adding accessories to the R82 Crocodile to either accommodate gait training with comprehensive support or walking with none or little support.

## Height adjustment

Easy, tool-free adjustment of the height of the top frame. This makes it easy to adjust the Crocodile, so it is a perfect fit for the user.

## Angle adjustable handles

The Crocodile has two adjustable handles. Positioning the handles close to the body gives better support and a better walking position.

## Wheel stop

Stopping the Crocodile from rolling backwards. Can easily be locked or unlocked.

## Folding mechanism

The frame of the Crocodile is foldable, making storage and transport easy.

## Tension adjustable rear wheels

Adjusting the tension on the rear wheels can make the resistance easier or harder for the individual user, according to their abilities.

## Front wheels with directional stabilizer

The possibility of locking the front wheels with the directional stabilisers, minimises the movement from side to side, making it easier for the child to control.



## Four sizes but five opportunities

The Crocodile is available in four sizes, but with the Starter Kit accessory, it is possible to retrofit the Crocodile size 1 to a smaller size, placing it between the Crocodile size 0 and size 1.




## Comparison for small users

Measurements		Size 0	Size 1 mounted with Starter Kit	Size 1
	Height from floor to handles	300 - 610 mm	380 - 700 mm	420 - 715 mm
	Width between handles	100 - 280 mm	90 - 350 mm	190 - 350 mm
	Length	640 mm	660 mm	720 mm



Size 0  
Available in white



Size 1  
Available in yellow

## Measurements

		Size 0	Size 1	Size 2	Size 3
	Width between handles	100-250 mm	190-350 mm	320-600 mm	330-640 mm
	Height from floor to handles	300-610 mm	420-715 mm 380-700 mm*	550-840 mm	750-1045 mm
	Max. user height	1050 mm	1150 mm	1500 mm	1800 mm
	Max. user mass	30 kg	30 kg	45 kg	80 kg
	Length	640 mm	720 mm	840 mm	870 mm
	Overall width	580 mm	620 mm	690 mm	690 mm
	Folded length	650 mm	720 mm	860 mm	1000 mm
	Folded width	580 mm	620 mm	690 mm	690 mm
	Folded height	300 mm	310 mm	270 mm	290 mm
	Turning width	1440 mm	1600 mm	1900 mm	2000 mm
	Mass device	6 kg	7 kg	7 kg	12 kg

\*Mounted with the Starter kit accessory. Further measurements on our website



Size 2  
Available in orange and metallic black



Size 3  
Available in anthracite grey

# Freedom of movement

The Crocodile encourages independent movement, exploration and discovery. It gives the child an opportunity to become more engaged and active within their environment and participate in activities with their peers and family.

The Crocodile gives children with limited movement the opportunity to be able to stand and walk. It allows freedom of movement and facilitates walking in the upright position.



## Support for the child

Additional support can be provided for posture and balance. This includes forearm support with handgrip, hip supports, back supports and side supports.

The Crocodile encourages the walking development for the user. It is versatile and adjustable and as the user develops or can walk with less support, the Crocodile allows you to easily reduce or adjust the support and challenge the child to achieve more.

The Crocodile challenges stepping, standing and balance development. It encourages children to move, become more independent and increases participation and engagement.



Follow me Dad. Lets go and  
feed the cute deer...



Watch our video with Konrad and his family having fun in the park. The Crocodile provides assistance for Konrad, so that he can move and explore the park independently in a safe and easy way.

# Meet Tobias and Oscar



## Tobias' story

When Tobias was about 9 months old, he was diagnosed with Spinal muscular atrophy type 1. Until a few years ago that would mean that children with this type of diagnosis, would experience that their muscles would become weak and waste away within a few years, making the child unable to meet typical milestones, like holding their head or sit independent. However, a new type of medicine has changed the life for children like Tobias.

After a few treatments Tobias could sit independently and a few months later he started using a R82 Pony for walking. Later he started using the R82 Mustang.



[Read the complete story about Tobias](#)

## Oscars story

Oscar has recently been diagnosed with unilateral right sided spastic cerebral palsy due to periventricular leukomalacia. Oscar lives in Denmark with his mother, father and older brother.

He has been classified at Gross Motor Function Classification System (GMFCS) level I or II and will be expected to be able to learn to walk independently between the age of 18 months to 2 years, without the need for any assistive devices.



[Read the complete story about Tobias](#)





# Meet Reid and Finn

## Reids story

Reid, 10 is a massive Derby County fan. He has played football since the age of three but left his local “able bodied” club after struggling with his pace. He began playing Frame Football at the Coventry club nearly two years ago and has since moved to newly established, Hucknall Rolls Royce, on home turf in Nottingham.

Reid’s Cerebral Palsy affects his mobility due to weakness in his legs. He is able to walk using sticks but uses a walker for more physical activities. As Reid explains “the Crocodile is the best frame by far! As it has a seat, I don’t have to take a wheelchair anywhere. I can walk and then take a break if I need to.



Watch Reid in action  
with his Crocodile



## Finns story

Finn, from Whitchurch, Bristol, began playing Frame Football in September 2015 and has progressed well in a short time.

Finn received a size 2 walker and said: “This is my first experience with a Crocodile Walker, it is the best frame I’ve ever used. I particularly like the seat as I can have a rest at half time. Usually, I play in goal but I have now tried playing in the outfield - awesome!”



See Finn playing in goal and  
outfield with his Crocodile

## Video guides


On our YouTube channel a wide range of videos are available - e.g. how to videos, product demonstrations and environmental videos.

Below is a selection of some of our how-to videos for the Crocodile.




 Surrounding support




 Grip bar on Crocodile 0




 Adjust the height and the handles



 Adjustments to the wheels



Visit our YouTube channel so see all of our Crocodile videos 



## Accessories




**Surrounding support**  
868610




**Gator back and side supports**  
86808, 86809

**Lateral strap**  
SN4016




**Back support**  
86886x, 86816




**Cross bar**  
86880



**Hip supports**  
86819



**Grip bar**  
86879-x



**Forearm support w/handgrip**  
95145-x



**Decelerating hand brake kit**  
86834-DEC, 86873-DEC

**Padded forearm positioner**  
869370

**Hand brake kit**  
86834, 86873




**Handle knobs**  
86861



**Flip down seat**  
8682x



**Sling seat**  
8682x



**Sling seat, light weight**  
8682x-1



**Anti-tips, fixed**  
86815-F, 86815-3F



**Attendant's handle**  
86866-1



**Starter kit**  
86801-0



**Bag**  
86839, 86843

**Anti-tips, swing-away**  
86815, 86815-3S

Compatibility	Size 0	Size 1	Size 2	Size 3
Surrounding support	•			
Gator back and side supports	•	•	•	•
Lateral strap		•	•	•
Back support	•	•	•	•
Cross bar				•
Hip supports		•	•	•
Grip bar	•	•	•	
Forearm support w/ handgrip	•	•	•	•
Padded forearm positioner	•	•	•	•
Decelerating hand brake kit		•	•	•
Hand brake kit		•	•	
Handle knobs	•	•	•	•
Flip down seat	•	•	•	•
Sling seat	•	•	•	
Sling seat light weight*	•	•	•	
Anti-tips, fixed	•	•	•	•
Anti-tips, swing-away	•	•	•	•
Attendant's handle	•	•	•	
Starter Kit		•		
Bag	•	•	•	•

\*Requires surrounding support or gator back and side support



Etac 2020 A/S. All rights reserved. Etac is a world-leading developer of innovative assistive devices and patient handling equipment. Since 1973, we have been committed to improving quality of life for the individual, family members and caregivers. For the intended purpose of the products and latest information, please visit [www.etac.com](http://www.etac.com). Etac is a registered trademark of Etac A/S. All rights reserved.

Etac is a world-leading developer of innovative assistive devices and patient handling equipment. Since 1973, we have been committed to improving quality of life for the individual, family members and caregivers.

For the intended purpose of the products and latest information, please visit [www.etac.com](http://www.etac.com)

R82 A/S

+45 79 68 58 88

[R82@etac.com](mailto:R82@etac.com)

[www.etac.com](http://www.etac.com)

**etac**<sup>®</sup>  
Creating Possibilities