


SOLUTIONS FOR THE FOOT AND ANKLE

Soft support edge with reduced pressure

SPORLASTIC AIR-MATRIX silicone pad: Anatomically shaped and nopped for pressure relief and proprioception

Innovative stretching zone for the perfect fit, even during movement



RELEVANT HERE:
Patient information for the ankle
Article No. 09310

Breathable 3D flat knitted fabric with recycled yarn

Three-part comfort zone for crease-free wearing

SUPREME series:
Physiotherapy to go
Access to therapy exercises and application instructions via QR code in the product packaging



MALLEO-Hi[®] SUPREME

PRODUCT AND INDICATION FINDER



	18	19	19	21	22	23	24	25	32	33
Page number										
Achillodynia (tendinopathy) with/without partial rupture							•	•		
Acute ankle injury				•		•				
Arthrosis	•	•	•		•					
Ligament sprain			•	•	•	•				
Bursitis subachillea, Bursitis calcanea							•	•		
Chronic capsular ligament instability	•	•	•	•	•	•				
Foot flexor paresis									•	•
Instabilities of the lower ankle			•	•	•	•				
Lumbar nerve root damage									•	•
Paratendinitis, paratendinopathy							•	•		
Prevention	•	•	•	•	•	•				
Quadriceps weakness										•
Achilles tendon reconstruction							•	•		
Rheumatism-related diseases of the ankle	•	•								
External / internal ligament rupture			•	•	•	•				
Swelling	•	•								
Ankle sprains	•	•	•	•	•	•				
Forefoot amputation (partial)										•

MALLEO-HIT® SUPREME

Ankle support with nopped AIR-MATRIX silicone friction pads



SPORLASTIC AIR-MATRIX silicone pad

Innovative stretching zone for the perfect fit, even during movement



Breathable 3D flat knitted fabric with recycled yarn

Three-part comfort zone for crease-free wearing



Physiotherapy to go
Access to therapy exercises and application instructions via QR code in the product packaging



MALLEO-HIT®

Ankle support with nopped silicone friction pads



friction pad

Non-elasticated lateral strap



friction pad

FIBULO-TAPE®

Functional ankle support as an alternative to tape bandaging



FOOT & ANKLE

ARTICLE NUMBER

07174

SIZE IN CM

18 Distal tibia circumference	
1	17 - 19
2	19 - 21
3	21 - 23
4	23 - 25
5	25 - 27
6	27 - 29

- Colours:
- Can be worn on the R and L

INDICATIONS

- Conservative / post-operative
- Chronic capsular ligament instability
- Mild ankle sprain
- Irritation of the upper / lower ankle
- Arthrosis of the upper / lower ankle
- Rheumatological conditions
- Swelling

EFFECT

- SPORLASTIC AIR-MATRIX silicone pad: Perforated, breathable, light, anatomically shaped and nopped silicone pad for pressure relief and proprioception
- Otherwise as with MALLEO-HIT®

TECHNIQUE

- Breathable 3D flat knitted fabric with recycled yarn
- Three-part comfort zone for crease-free wearing
- Soft support edge with reduced pressure
- Innovative stretching zone for the perfect fit, even during movement

ARTICLE NUMBER

07074

SIZE IN CM

18 Distal tibia circumference	
1	17 - 19
2	19 - 21
3	21 - 23
4	23 - 25
5	25 - 27
6	27 - 29

- Colours:
- Can be worn on the R and L

KIDS

87074

18 Distal tibia circumference	
C	16 - 18

- Colour:
- Can be worn on the R and L

INDICATIONS

- Conservative / post-operative
- Chronic capsular ligament instability
- Mild ankle sprain
- Irritation of the upper / lower ankle
- Arthrosis of the upper / lower ankle
- Rheumatological conditions
- Swelling

EFFECT/TECHNIQUE

- Stabilisation, compression and relief
- Optimisation of proprioception
- Pain relief through reduced nociceptive afferences
- Friction massage
- Accelerated reabsorption of effusions
- Beneficial effects on insertion tendinopathies and bursitis typically seen with arthrosis
- Extremely elastic specialist knit for virtually crease-free wearing

ARTICLE NUMBER

07073

SIZE IN CM

18 Distal tibia circumference	
1	17 - 19
2	19 - 21
3	21 - 23
4	23 - 25
5	25 - 27
6	27 - 29

- Colours:
- Please specify R or L

INDICATIONS

- Conservative / post-operative
- Acute capsular ligament rupture
- Chronic capsular ligament instability
- Primary and secondary prophylaxis of sprains, especially during sport
- Ankle arthrosis
- Stabilisation of the upper ankle

EFFECT/TECHNIQUE

- Protective effect like tape bandaging
- Stabilisation of the anterior talo-fibular ligament and fibulo-calcaneal ligament
- Reduction of the risk of supination trauma
- Optimisation of proprioception
- Gentle on the skin
- 8 straps
- Non-elasticated lateral strap