

MANU-HiT® POLLEX

Wrist brace with thumb enclosure



MANU-HiT®

Wrist brace with volar and dorsal reinforcement



Mouldable aluminium inlay



Custom-mouldable volar and dorsal reinforcement bars

ARTICLE NUMBER

07645

SIZE IN CM

6 Wrist circumference

XS	13 - 15
S	15 - 17
M	17 - 19
L	19 - 22
XL	22 - 25

- Length: 18 cm
- Colours:
- Please specify R or L

INDICATIONS

- Conservative / post-operative
- Injuries to the carpal bones
- Arthrosis of the radial parts of the carpus
- Arthrosis, sprain, ulnar lateral ligament lesion of the metacarpo-phalangeal joint
- Scaphoid injuries
- Carpal tunnel syndrome
- Radial instability of the wrist
- Peripheral radial nerve paralysis with drop hand and poor abduction of the thumb

EFFECT/TECHNIQUE

- Stabilisation and partial immobilisation through thumb enclosure in a functional position
- Preservation of grip function
- Reduction of strain on the carpal bones, joint capsule and ligaments
- Anatomically pre-moulded reinforcement inlay
- SOFT-EDGE TECHNOLOGY: Soft material blend for maximum wearing comfort

ARTICLE NUMBER

07035

SIZE IN CM

6 Wrist circumference

XS	13 - 15
S	15 - 17
M	17 - 19
L	19 - 22
XL	22 - 25

- Length: 18 cm
- Colours:
- Please specify R or L

KIDS

87035

6 Wrist circumference

C	12 - 14
---	---------

- Length: 17 cm
- Colour:
- Please specify R or L

INDICATIONS

- Conservative / post-operative
- Following sprains
- Tendon sheath inflammation
- Arthrosis of the wrist
- Wrist strain symptoms
- Injuries to the carpal bones
- Following removal of plaster casts
- Navicular pseudo-arthritis
- Drop hand
- Lunatomalacia
- Rheumatoid arthritis
- Carpal tunnel syndrome

EFFECT/TECHNIQUE

- Reinforced stabilisation and restriction of movement
- Reduction of strain on the carpal bones, joint capsule and ligaments
- Closely fits the body's shape
- Fixation of the carpus, free mobility of the thumb and fingers
- Custom-mouldable volar and dorsal reinforcement bars
- SOFT-EDGE TECHNOLOGY: Soft material blend for maximum wearing comfort

MANU-HiT® CARPAL

Carpal tunnel syndrome brace with CTS RELIEF PAD to protect against pressure irritation



SOFT-EDGE TECHNOLOGY:
Soft material blend for maximum wearing comfort



CTS RELIEF PAD:
Anatomically shaped with exposure of the carpal tunnel and the median nerve



Physiotherapy to go
Access to therapy exercises and application instructions via QR code in the product packaging



ARTICLE NUMBER

07233

SIZE IN CM

6 Wrist circumference

XS	13 - 15
S	15 - 17
M	17 - 19
L	19 - 22
XL	22 - 25

- Length: 18 cm
- Colours:
- Please specify R or L

INDICATIONS

- Conservative / post-operative
- Carpal tunnel syndrome/ CTS
- Tendon sheath inflammation
- Arthrosis of the wrist
- Peripheral irritation of the ulnar nerve in the Guyon canal
- Post-operative immobilisation
- Post-operative residual symptoms

EFFECT

- Relief of carpal tunnel
- Fixation in functional position
- Avoidance of further narrowing of the carpal tunnel cross-section
- Protection against pressure irritation in the area of the carpal tunnel
- Stimulation of oedema drainage in the wrist and vascular system

TECHNIQUE

- CTS RELIEF PAD: anatomically shaped with exposure of the carpal tunnel, the median nerve, one-piece, individually positionable via micro-Velcro and removable
- Closure system: Three supple micro-Velcro fasteners and slip stop
- Optimum anatomical fit, also possible with conical forearm shapes
- Custom-adjustable thumb fastener
- Custom-mouldable volar and dorsal reinforcement bars

MANU-HiT® CLASSIC

Wrist brace with volar and dorsal reinforcement



MANU-HiT® POLLEX CLASSIC

Wrist brace with thumb enclosure



Custom-mouldable volar and dorsal reinforcement bars



ARTICLE NUMBER

07031

Length: 18 cm

07032

Length: 22 cm

SIZE IN CM

6 Wrist circumference

XS	13 - 15
S	15 - 17
M	17 - 19
L	19 - 22
XL	22 - 24

XS only with Art. No. 07031

INDICATIONS

- Conservative / post-operative
- Wrist sprain trauma
- Tendon sheath inflammation
- Arthrosis of the wrist
- Strain symptoms
- Following removal of plaster casts
- Navicular pseudo-arthritis
- Drop hand
- Lunatomalacia
- Rheumatoid arthritis
- Carpal tunnel syndrome

EFFECT/TECHNIQUE

- Increased stabilisation and restriction of movement
- Reduction of strain on the carpal bones, joint capsule and ligaments
- Closely fits the body's shape
- Fixation of the carpus, free mobility of the thumbs and fingers
- **Article No. 07031**
- Colours:
- Please specify R or L
- **Article No. 07032**
- Colour:
- Please specify R or L

ARTICLE NUMBER

07640

INDICATIONS

- Conservative / post-operative
- Arthrosis of the radial parts of the carpus
- Arthrosis, sprain, ulnar lateral ligament lesion of the metacarpo-phalangeal joint
- Scaphoid injuries
- Carpal tunnel syndrome
- Radial instability of the wrist
- Peripheral radial nerve paralysis with drop hand and poor abduction of the thumb

EFFECT/TECHNIQUE

- Stabilisation and partial immobilisation through thumb enclosure in a functional position
- Preservation of grip function
- Reduction of strain on the carpal bones, joint capsule and ligaments
- Anatomically pre-moulded reinforcement inlay